

Product datasheet

Recombinant human DGKA protein (Active) ab268441

2 Images

Description

Product name	Recombinant human DGKA protein (Active)
Biological activity	The specific activity of ab268441 was 29.1 nmol/min/mg in a kinase assay using Dilauroyl-sn-glycerol as substrate.
Purity	> 90 % SDS-PAGE.
Expression system	Baculovirus infected Sf9 cells
Accession	P23743
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>MAKERGLISPSDFAQLQKYMESTKKVSDVLKLFEDGEM AKYVQGDAIGY EGFQQFLKMLEVDNVPRHLSLALFQSFETGHCLNETNVT KDVVCLNDVS CYFSLLEGGRPEDKLEFTFKLYDTRNGILDSSEVDKILQ MMRVAEYLD WDVSELRPILQEMMKEIDYDGGSGSVSQA EWVRAGATTVP LLVLLGLEMTL KDDGQHMWRPKRFRPVYCNLCESSIGLGKQGLSCNLC KYTVHDQCAMKA LPCEVSTYAKSRKDIGVQSHVWRGGCESGRCDRCQKKI RYHSLTGLHC VWCHLEIHDDCLQAVGHECDCGLLRDHILPPSSMPSVLA SGPDRKNSKT SQKTMDLNLSTSEALRIDPVPNTHPLLVFVNPKSGGKQG QRVLWKFQYI LNPRQVFNLLKDGPEIGLRLFKDVPDSRILVCGGDGTVGW ILETIDKANL PVLPPVAVLPLGTGNDLARCLRWGGGYEQNLAKILKDL EMSKVHMDRW SVEVIPQQTEEKSDPVFQIINNYFSIGVDASIAHRFHIMRE KYPEKFNS RMKNKLWYFEFATSESIFFSTCKKLEESLVEICGKPLDLSN</p>

LSLEGI AVL
NIPSMHGGSNLWGDTRRPHGDIYGINQALGATAKVITDPDIL
KTCVPDLS
DKRLEVVGLEGAIEMGQIYTKLKNAGRRLAKCSEITFHHTK
TLPMQIDGE
PWMQTPCTKITHKNQMPMLMGPPPRSTNFFGFLS

Molecular weight information	~82 kDa by SDS-PAGE
Amino acids	1 to 735
Tags	His tag N-Terminus
Additional sequence information	GenBank: NM_001345

Specifications

Our [Abpromise guarantee](#) covers the use of **ab268441** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Functional Studies SDS-PAGE
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. pH: 7.00 Preservative: 1.02% Imidazole Constituents: 0.82% Sodium phosphate, 1.74% Sodium chloride, 0.002% PMSF, 0.004% DTT, 25% Glycerol (glycerin, glycerine) This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------------------------	--

General Info


Function	Upon cell stimulation converts the second messenger diacylglycerol into phosphatidate, initiating the resynthesis of phosphatidylinositols and attenuating protein kinase C activity.
Tissue specificity	Lymphocytes and oligodendroglial cells.
Sequence similarities	Belongs to the eukaryotic diacylglycerol kinase family. Contains 1 DAGKc domain. Contains 2 EF-hand domains. Contains 2 phorbol-ester/DAG-type zinc fingers.

Images



The specific activity of ab268441 was 29.1 nmol/min/mg in a kinase assay using Dilauroyl-sn-glycerol as substrate.

Functional Studies - Recombinant human DGKA protein (Active) (ab268441)



SDS-PAGE analysis of ab268441.

SDS-PAGE - Recombinant human DGKA protein (Active) (ab268441)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors