

Product datasheet

Recombinant Human DPF2/REQ protein ab159327

1 Image

Description

<b>Product name</b>	Recombinant Human DPF2/REQ protein
<b>Expression system</b>	Wheat germ
<b>Protein length</b>	Full length protein
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Sequence</b>	<p>MAAVVENVVKLLGEQYK DAMEQCHNYNARLCAERSVRL  PFLDSQTGVAQ  SNCYWMEKRHRGPGGLASGQLYSYPARRWRKKRAHPP  EDPRLSFPSIKP  DTDQTLKKEGLISQDGSSLEALLRTDPLEKRGAPDPRVD  DDSLGEFPVTN  SRARKRILEPDDFLDDLDEDEYEDTPKRRGKKGSKGKG  VGSARKKLDAS  ILED RDKPYACDICGKRYKNRPGLSYHYAHSHLAEEEGED  KEDSQPPTPV  SQRSEEQSKKGPDGLALPNNYCDFCLGDSKINKKTGQP  EELVSCSDCGR  SGHP SCLQFTPVMMAAVKTYRWQCIECKCCNICGTSEND  DQLLFCDDCDR  GYHMYCLTPSMSEPPEGSWSCHLCLDLLKEKASYQNQN  SS</p>
<b>Amino acids</b>	1 to 391
<b>Tags</b>	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab159327** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	ELISA Western blot
<b>Form</b>	Liquid

## Additional notes

This product was previously labelled as DPF2.

## Preparation and Storage

### Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

## General Info

### Function

May be a transcription factor required for the apoptosis response following survival factor withdrawal from myeloid cells. Might also have a role in the development and maturation of lymphoid cells.

### Tissue specificity

Ubiquitous.

### Sequence similarities

Belongs to the requiem/DPF family.


Contains 1 C2H2-type zinc finger.

Contains 2 PHD-type zinc fingers.

### Cellular localization

Nucleus. Cytoplasm. 30% nuclear. 70% cytoplasmic.

## Images



ab159327 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human DPF2/REQ protein (ab159327)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors