

Product datasheet

Recombinant Human DPH3P1 protein ab177714

1 Image

Description

Product name	Recombinant Human DPH3P1 protein	
Purity	> 90 % SDS-PAGE. ab177714 was purified using conventional chromatography.	
Expression system	Escherichia coli	
Accession	Q9H4G8	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMVVFHDE VEIEDFQYDE DSETYFCPCP CGDNFSITKE ELENAGEVAM CPGCSLIKV YDKDQFACG ETVPVPSVVK E	
Predicted molecular weight	11 kDa including tags	
Amino acids	1 to 78	
Tags	His tag N-Terminus	
Additional sequence information	NP_542788	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab177714** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------------------------	--

pH: 8.00


Constituents: 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine)

General Info

Relevance

Belongs to the DPH3 family.

Images



15% SDS-PAGE analysis of ab177714 (3µg).

SDS-PAGE - Recombinant Human DPH3P1 protein
(ab177714)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors