

Product datasheet

Recombinant Human DUS1L protein ab163785

1 Image

Description

Product name	Recombinant Human DUS1L protein	
Expression system	Wheat germ	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MPKLQGFEFWSRTLGRARHVVAPMVDQSELAWRLLS RRHGAQLCYTPMLH AQVFVRDANYRKENLYCEVCPEDRPLVQFCANDPEV FVQAALLAQDYCD AIDLNLGCPQMIAKRGHYGAFLQDEWDLLQRMILLAHE KLSVPVTCKIRV FPEIDKTVRYAQMLEKAGCQLLTVHGRTKEQKGPLSG AASWEHIKAVRKA VAIPVFANGNIQCLQDVERCLRDTGVQGVMSAEGNLH NPALFEGRSPAVW ELAEEYLDIVREHPCPLSYVRAHLFKLWHHTLQVHQEL REELAKVKTLEG IAVSQELKLRQCQEEISRQEGAKPTGDLPFHWICQPYR PGPREGSKEKA GARSKRALEEEEGGTEVLSKNKQKKQLRNPHKTFDP SLKPKYAKCDQCGN PKGNRCVFSLCRGCCCKRASKETADCPGHLLFKTK LEKSLAWKEAQPEL QEPQPAAPGTPGGFSEVMGSALA	
Amino acids	1 to 473	
Tags	GST tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab163785** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
---------------------	-------

Western blot

Form Liquid

Additional notes

Preparation and Storage


Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function Catalyzes the synthesis of dihydrouridine, a modified base found in the D-loop of most tRNAs.

Sequence similarities Belongs to the Dus family. Dus1 subfamily.

Images



ab163785 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human DUS1L protein
(ab163785)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors