

Product datasheet

Recombinant Human ECHDC1 protein ab123187

1 Image

Overview

Product name	Recombinant Human ECHDC1 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Escherichia coli
Amino Acid Sequence	
Accession	Q9NTX5
Species	Human

Sequence	<p>MGSSHHHHHH SGLVPRGSH MALKQEMAKS LLKTASLSGR TKLLHQTGLS LYSTSHGFYE EEVKKTLQQF PGGSIDLQKE DNGIGILTLN NPSRMNAFSG VMMLQLEKV IELENWTEGK GLIVRGAKNT FSSGSDLNAV KSLGTPEDGM AVCMFMQNTL TRFMRLPLIS VALVQGWALG GGAEFTTACD FRLMTPESKI RFVHKEMGII PSWGGTTRLV EIIGSRQALK VLSGALKLDS KNALNIGMVE EVLQSSDETK SLEEAQEWLK QFIQGPPEVI RALKKSVCSG RELYLEEALQ NERDLLGTVW GGPANLEAIA KKGKFNK</p>
-----------------	--

Molecular weight	36 kDa including tags
Amino acids	1 to 307
Tags	His tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab123187** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	<p>SDS-PAGE</p> <p>Mass Spectrometry</p>
Mass spectrometry	MALDI-TOF

Purity > 95 % SDS-PAGE.
ab123187 is purified using conventional chromatography techniques.

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol, 0.58% Sodium chloride

General Info

Sequence similarities Belongs to the enoyl-CoA hydratase/isomerase family.

Images



15% SDS-PAGE showing ab123187 (3µg).

SDS-PAGE - Recombinant Human ECHDC1 protein
(ab123187)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors