

Product datasheet

Recombinant Human ERGIC3 protein ab127471

Overview

Product name	Recombinant Human ERGIC3 protein
Protein length	Protein fragment

Description

Nature	Recombinant
Source	Escherichia coli
Amino Acid Sequence	
Accession	Q9Y282
Species	Human
Molecular weight	30 kDa
Amino acids	78 to 343
Tags	His-DHFR tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab127471** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Purity	= 90 % Purified via His tag. Purity determined to be 90% by SDS-PAGE.
Form	Lyophilised

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
Reconstitution	Reconstitute with water to desired concentration

General Info

Function	Possible role in transport between endoplasmic reticulum and Golgi.
Sequence similarities	Belongs to the ERGIC family.
Cellular localization	Endoplasmic reticulum-Golgi intermediate compartment membrane. Golgi apparatus > cis-Golgi network membrane. Endoplasmic reticulum membrane. Cycles between the endoplasmic reticulum and the Golgi.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors