

Product datasheet

Recombinant human FAK protein ab105909

[1 References](#) [2 Images](#)

Overview

Product name	Recombinant human FAK protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Baculovirus infected Sf9 cells
Amino Acid Sequence	
Accession	Q05397
Species	Human
Molecular weight	140 kDa
Amino acids	1 to 1052
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab105909** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Biological activity	The Specific activity of ab105909 was determined to be 28 nmol/min/mg.
Applications	Functional Studies SDS-PAGE
Purity	> 85 % SDS-PAGE. Purity was determined to be >85% by densitometry.
Form	Liquid
Additional notes	ab204877 (Poly (4:1 Glu, Tyr) peptide) can be utilized as a substrate for assessing kinase activity

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.307% Glutathione, 0.01% Triton-X-100, 0.00174% PMSF, 0.00385% DTT, 0.79%
------------------------------	---

Tris HCl, 0.00292% EDTA, 25% Glycerol, 0.87% Sodium chloride

This product is an active protein and may elicit a biological response in vivo, handle with caution.

General Info

Function

Non-receptor protein-tyrosine kinase implicated in signaling pathways involved in cell motility, proliferation and apoptosis. Activated by tyrosine-phosphorylation in response to either integrin clustering induced by cell adhesion or antibody cross-linking, or via G-protein coupled receptor (GPCR) occupancy by ligands such as bombesin or lysophosphatidic acid, or via LDL receptor occupancy. Microtubule-induced dephosphorylation at Tyr-397 is crucial for the induction of focal adhesion disassembly. Plays a potential role in oncogenic transformations resulting in increased kinase activity.

Tissue specificity

Expressed in all organs tested, in lymphoid cell lines, but most abundantly in brain.

Sequence similarities

Belongs to the protein kinase superfamily. Tyr protein kinase family. FAK subfamily.

Contains 1 FERM domain.

Contains 1 protein kinase domain.

Domain

The first Pro-rich domain interacts with the SH3 domain of CRK-associated substrate (BCAR1) and CASL.

The carboxy-terminal region is the site of focal adhesion targeting (FAT) sequence which mediates the localization of FAK1 to focal adhesions.


Post-translational modifications

Phosphorylated on 6 tyrosine residues upon activation. Microtubule-induced dephosphorylation at Tyr-397 could be catalyzed by PTPN11 and regulated by ZFYVE21. Dephosphorylated by PTPN11 upon EPHA2 activation by its ligand EFNA1.

Cellular localization


Cell junction > focal adhesion. Cell membrane. Constituent of focal adhesions.

Images



Kinase Assay demonstrating specific activity of ab105909.

Functional Studies - Recombinant human FAK protein (ab105909)



SDS-PAGE showing ab105909 at approximately 140kDa.

SDS-PAGE - Recombinant human FAK protein (ab105909)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors