

Product datasheet

Recombinant Human FLJ36353 protein (His tag) ab196051

Description

Product name	Recombinant Human FLJ36353 protein (His tag)	
Purity	> 95 % SDS-PAGE. SEC-HPLC and reducing SDS-PAGE	
Endotoxin level	< 1.000 Eu/μg	
Expression system	HEK 293 cells	
Accession	Q9BSG5	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	SRPLQARSQQHHGLAADLGKGLHLAGPCCPSEMDTTE TSGPGNHPERCG VPSPECESFLEHLQR ALRSRFRLLRLLGVRQAQPLCEELCQAWFANCEDDITCGP TWLPLSEKRGCEPSCLTYGQTFADGT DLCRSALGHALPVAAPGARHCFNISISAVPRPRPGRRGRE APSRRRSRSPR TSILDAAGSGSGSGS GSGPVDHHHHHH*	
Predicted molecular weight	22 kDa including tags	
Amino acids	31 to 229	
Tags	His tag C-Terminus	
Additional sequence information	This product is for the mature full length protein from aa 31 to 229. The signal peptide is not included.	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab196051** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	HPLC SDS-PAGE
Form	Lyophilized

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at -20°C.

pH: 7.4

Constituents: 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride, 0.02% DTT

Lyophilized from a 0.2 µm filtered solution

General Info

Tissue specificity

Expressed in retina.

Sequence similarities

Belongs to the folate receptor family.

Cellular localization

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors