

Product datasheet

Recombinant Human FLNA protein ab131760

1 Image

Description

Product name	Recombinant Human FLNA protein
Expression system	Wheat germ
Accession	Q96C61
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human

Sequence

```

MPSGKVAQPTITDNKDGT/TVRYAPSEAGLHEMDIRYDNM
HIPGSPLQFY
VDYVNCGHVTA YGPGLTHGVVNKPATFTVNTKDAGEGGL
SLAIEGPSKAE
ISCTDNQDGTCSVSYLPVLPGDYSILVKYNEQHVP GSPFT
ARVTGDDSMR
MSHLKVGSAADIPINSETDLSLLTATVPPSGREEPCLLK
RLRNGHVG I
SFVPKETGEHLVHVKKNQHVASSPIPVVISQSEIGDASR
VRVSGQQLHE
GHTFEPAEFIIDTRDAGYGGLSLSIEGPSKVDINTEDLEDGT
CRVTYCPT
EPGNYIINIKFADQHVP GSPFSVKVTGEGRVKESITRRRRA
PSVANVGSH
CDLSLKIPEISIQDMTAQVTSPSGKTHEAEIVEGENHTYCIR
FVPAEMGT
HTVSVKYGQHVPGSPFQFTVGPLGEGGAHKVRAGGPG
LERA EAGVPAEF
SMTREAGAGGLAIAVEGPSKAEISFEDRKDGSCGVAYV
QEPGDYEVS
KFNEEHIPDSPFVVPVASPSGDARRLTVSSLQESGLKVN
QPASFAVSLNG
AKGAIDAKVHSPSGALEECYVTEIDQDKYAVRFIPRENGVY
LIDVKFNGT
HIPGSPFKIRVGE PGHGGDPGLVSAYGAGLEGGVTGNPA
EFVVNTSNAGA
GALSVTIDGPSKVKMDCQECPEGYRVTYTPMAPGSYLISIK
YGGPYHIGG

```

SPFKAKVTGPRLVSNHSLHETSSVFVDSLTKATCAPQHG
APGPGPADASK
VVAKGLGLSKAYVQKSSFTVDCSKAGNNMLLVGVHGP
RTPCEEILVKHV
GSRLYSVSYLLKDKGEYTLVVKWGDEHIPGSPYRVVVP

Predicted molecular weight	118 kDa including tags
Amino acids	1 to 838
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab131760** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot SDS-PAGE
---------------------	-----------------------------------


Form	Liquid
-------------	--------

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Images



12.5% SDS-PAGE stained with Coomassie Blue showing ab131760.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors