

Product datasheet

Recombinant Human FRMD8 protein ab132735

1 Image

Description

Product name	Recombinant Human FRMD8 protein
Expression system	Wheat germ
Accession	Q9BZ67
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>MDGTEGSAGQPGAERSHRSSVSSVGARAADVLYL ADDTVVPLAVENLP SLSAHELHRAVREVLQLPDIALDVFALWLVSPLEVLQ KPKHQPYKLGRQ WPELLLRFTSAPDDDAMDEPFLQFRRNVFFPKRRE LQIHDEEVLRLLYE EAKGNVLAARYPCDVEDCEALGALVCRVQLGPYQPG RPAACDLREKLDSF LPAHLCKRGQSLFAALRGRGARAGPGEQGLLNAYRQ VQEVSSDGGCEAAL GTHYRAYLLKCHELPHYGCAFFHGEVDKPAQGFLHRG GRKPVSV AISLEG VHVIDSREKHVLLGLRFQELSWDHTSPEEEEEPIWLEF DGDSEGTPVKNL LKIYSKQAEMLMSSLIEYCIELSQA AEPAGPQDSATGSPS DPSSSLAPVQR PKLRRQGSVSSRIQHLSTIDYVEDGKGIRRVKPKRTTS FFSRQLSLGQG SYTVVQPGDSLEQG</p>
Predicted molecular weight	78 kDa including tags
Amino acids	1 to 464

Specifications

Our [Abpromise guarantee](#) covers the use of **ab132735** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Western blot

	SDS-PAGE
	ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 µg/µl.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------

General Info

Sequence similarities	Contains 1 FERM domain.
------------------------------	-------------------------

Images



12.5% SDS-PAGE analysis of ab132735 stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors