

Product datasheet

Recombinant Human FSBP protein (denatured) ab156713

1 Image

Overview

| | |
|-----------------------|--|
| Product name | Recombinant Human FSBP protein (denatured) |
| Protein length | Full length protein |
| Description | Recombinant Human FSBP protein |

Description

| | |
|----------------------------|---|
| Nature | Recombinant |
| Source | Escherichia coli |
| Amino Acid Sequence | |
| Accession | O95073 |
| Species | Human |
| Sequence | <p>MGSSHHHHHHSSGLVPRGSHMGSMVGKARSSNFTLS EKDLLKLVKPYVK ILEEHTNKHSVIVEKNRCWDIIAVNYNAIGVDRPPRTAQ GLRTLKRLKE YAKQELLQQKETQSDFKSNISEPTKKVMEMIPQISSFC LVRDRNHIQSAN LDEEAQAGTSSLQVMLDHPVAITVEVKQEEDIKPPP PLVLNSQQSDTLE QREEHELHVHMERSLSPLSSVDMRMTSSPSSIPRRD DFFRHESGEHFRS LLGYDPQILQMLKEEHQILENQKNFGLYWQEKRDLKR RQQLLEEELLRA KIEVEKLRKAIRLRHDLPYNSL</p> |
| Molecular weight | 37 kDa including tags |
| Amino acids | 1 to 299 |
| Tags | His tag N-Terminus |

Specifications

Our [Abpromise guarantee](#) covers the use of **ab156713** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|---------------------|------------------|
| Applications | SDS-PAGE |
| Purity | > 85 % SDS-PAGE. |
| Form | Liquid |


Preparation and Storage

| | |
|------------------------------|---|
| Stability and Storage | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 6% Urea, 0.32% Tris HCl, 10% Glycerol |
|------------------------------|---|

General Info

| | |
|------------------------------|---------|
| Cellular localization | Nuclear |
|------------------------------|---------|

Images



15% SDS-PAGE analysis of ab156713 (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors