

Product datasheet

Recombinant Human FUNDC1 protein ab165435

1 Image

Description

Product name	Recombinant Human FUNDC1 protein	
Expression system	Wheat germ	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MATRNPPPQDYESDDDSYEVLDLTEYARRHQWWNRVFG HSSGPMVEKYSV ATQVMGGVTGWCAGFLFQKVGKLAATAVGGGFLLLQIAS HSGYVQIDWK RVEKDVNKAKRQIKKRANKAAPEINNLEEATEFIKQNVISS GFVGGFL LGLAS	
Amino acids	1 to 155	
Tags	GST tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab165435** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot
Form	Liquid

Additional notes


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
	pH: 8.00
	Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function	Acts as an activator of hypoxia-induced mitophagy, an important mechanism for mitochondrial quality control.
Tissue specificity	Widely expressed.
Sequence similarities	Belongs to the FUN14 family.
Domain	The YXXL motif mediates the interaction with MAP1 LC3 family proteins MAP1LC3A, MAP1LC3B and GABARAP (PubMed:22267086).
Post-translational modifications	Phosphorylation at Tyr-18 by SRC inhibits activation of mitophagy. Following hypoxia, dephosphorylated at Tyr-18, leading to interaction with MAP1 LC3 family proteins and triggering mitophagy.
Cellular localization	Mitochondrion outer membrane.

Images



ab165435 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors