

Product datasheet

Recombinant Human Galectin 2 protein ab99239

1 Image

Description

Product name	Recombinant Human Galectin 2 protein	
Purity	> 95 % SDS-PAGE. ab99239 is purified using conventional chromatography techniques.	
Expression system	Escherichia coli	
Accession	P05162	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MGSSHHHHHSSGLVPRGSHMTGELEVKNMDMKPGST LKITGSIADGTDG FVINLGQGTDKLNLHFNPRFSESTVCMNSLDGSNWGQEQR EDHLCFSPGS EVKFTVTFESDKFKVKLPDGHELTFPNRLGSHLSYLSVR GGFNMSSFKL KE	
Predicted molecular weight	17 kDa including tags	
Amino acids	1 to 132	
Tags	His tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab99239** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw

cycles.

pH: 8.00


Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 1.16% Sodium chloride

General Info

Function This protein binds beta-galactoside. Its physiological function is not yet known.

Sequence similarities Contains 1 galectin domain.

Images



15% SDS-PAGE showing ab99239 at approximately 16.8kDa (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors