

Product datasheet

Recombinant Human HOXD3 protein ab158688

1 Image

Overview

Product name	Recombinant Human HOXD3 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human
Sequence	<p>MLFEQQQQAELPECTMQKAAYENPGLFGGYGYSKT TDTYGYSTPHQPY PPPAAASSLDTDYPGSACSIQSSAPLRAPAHKGAELN GSCMRPGTGNSQG GGGGSQPPGLNSEQQPPQPPPPPTLPPSSPTNPGG GVPAKKPKGGPNAS SSSATISKQIFPWMKESRQNSKQKNSCATAGESCEDK SPPGPASKRVRTA YTSAQLVELEKEFHFNRYLCRPRRVEMANLLNLTERQI KIWFQNRMMKYK KDQKAKGILHSPASQSPERSPPLGGAAGHVAYSQQLP PVPGLAYDAPSPP AFAKSQPNMYGLAAYTAPLSSCLPQQKRYAAPEFEPH PMASNGGGFASAN LQGSPVYVGGNFVESMAPASGPVFNLGHLSHPSAS VDYSCAAQIPGNHH HGPCDPHPTYTDLSAHHSSQGRLPEAPKLTHL</p>
Amino acids	1 to 432
Tags	proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab158688** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
---------------------	-------

	Western blot
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Function	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
Sequence similarities	Belongs to the Antp homeobox family. Contains 1 homeobox DNA-binding domain.
Cellular localization	Nucleus.

Images



ab158688 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human HOXD3 protein
(ab158688)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors