

Product datasheet

Recombinant Human HSCB protein ab128453

1 Image

Description

Product name	Recombinant Human HSCB protein
Purity	> 90 % SDS-PAGE. ab128453 was purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	Q8MWL3
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSHEMAASQA GSNYPRCWNC GGPWGPGRED RFFCPQCRAL QAPDPTRDYF SLMDCNRSFR VDTAKLQHRY QQLQRLVHPD FFSQRSQTEK DFSEKHSTLV NDAYKTLLAP LSRGLYLLKL HGIEIPERTD YEMDRQFLIE IMEINEKLAE AESEAAMKEI ESIVKAKQKE FTDNVSSAFE QDDFEEAKEI LTKMRYFSNI EEKIKLKKIP L
Predicted molecular weight	27 kDa including tags
Amino acids	30 to 235
Tags	His tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab128453** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.88% Sodium chloride

General Info

Function

May act as a co-chaperone in iron-sulfur cluster assembly in mitochondria.

Tissue specificity

Liver, muscle and heart.

Pathway

Cofactor biosynthesis; iron-sulfur cluster biosynthesis.

Sequence similarities


Belongs to the hscB family.

Contains 1 J domain.

Cellular localization

Mitochondrion.

Images



15% SDS-PAGE analysis of 3 µg ab128453.

SDS-PAGE - Recombinant Human HSCB protein
(ab128453)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors