

Product datasheet

Recombinant Human Hsp47 protein ab86918

1 Image

Description

<b>Product name</b>	Recombinant Human Hsp47 protein
<b>Purity</b>	> 90 % SDS-PAGE. ab86918 was purified using conventional chromatography techniques.
<b>Expression system</b>	Escherichia coli
<b>Protein length</b>	Full length protein
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Sequence</b>	<p><b>M</b>RGSHHHHHH <b>G</b>MASMTGGQQ <b>M</b>GRDLYDDDD  <b>K</b>DRWGSMAA EVKKPAAAAA PGTAEKLSPK  AATLAERSAG LAFSLYQAMA KDQAVENILV  SPVVVASSLG LVSLGGKATT ASQAKAVLSA  EQLRDEEVHA GLGELLRSLN NSTARNVTWK  LGSRLYGPSS VSFADDFVRS SKQHYNCEHS  KINFRDKRSA LQSINEWAAQ TTDGKLPEVT  KDVERTDGAL LVNAMFFKPH WDEKFHHKMV  DNRGFMVTRS YTVGVMMHR TGLYNYDDE  KEKLQVEMP LAHKLSSLII LMPHHVEPLE RLEKLLTKEQ  LKWMGKMQK KAVAIPLPKG VVEVTHDLQK  HLAGLGLTEA IDKNKADLSR MSGKKDLYLA  SVFHATAFEL DTDGNPFDQD IYGREELRSP KLFYADHPFI  FLVRDTQSGS LLFIGRLVRP KGDKMRDEL</p>

Specifications

Our [Abpromise guarantee](#) covers the use of **ab86918** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** SDS-PAGE

**Form** Liquid

Preparation and Storage

## Stability and Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

## General Info

---

### Function

Binds specifically to collagen. Could be involved as a chaperone in the biosynthetic pathway of collagen.

### Involvement in disease

Note=Defects in SERPINH1 may cause severe autosomal recessive osteogenesis imperfecta (OI). Osteogenesis imperfecta defines a group of connective tissue disorders characterized by bone fragility and low bone mass.

### Sequence similarities

Belongs to the serpin family.

### Cellular localization

Endoplasmic reticulum lumen.

---

## Images

---



15% SDS-PAGE analysis of 3µg ab86918

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors