

Product datasheet

Recombinant human HSPA2 protein ab92433

2 Images

Description

Product name	Recombinant human HSPA2 protein
Biological activity	ATPase activity: Positive
Purity	> 95 % SDS-PAGE. Purified by Multi-step chromatography. > 95% pure as determined by SDS-PAGE and Western Blot analysis.
Expression system	Escherichia coli
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human

Specifications

Our [Abpromise guarantee](#) covers the use of **ab92433** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot Functional Studies SDS-PAGE
Form	Liquid
Additional notes	Endotoxin: >500 EU/mg as determined by LAL gel clot assay

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------------------------	---

General Info


Function	In cooperation with other chaperones, Hsp70s stabilize preexistent proteins against aggregation and mediate the folding of newly translated polypeptides in the cytosol as well as within
-----------------	---

organelles. These chaperones participate in all these processes through their ability to recognize nonnative conformations of other proteins. They bind extended peptide segments with a net hydrophobic character exposed by polypeptides during translation and membrane translocation, or following stress-induced damage.

Sequence similarities

Belongs to the heat shock protein 70 family.

Images



SDS-PAGE - Recombinant human HSPA2 protein (ab92433)


SDS-PAGE showing ab92433:

Lane 1: 0.5µg

Lane 2: 1µg

Lane 3: 2µg

Lane 4: 5µg



Western blot - Recombinant human HSPA2 protein (ab92433)

All lanes : anti-HSPA2

Lane 1 : Recombinant human HSPA2 protein (ab92433)

Lane 2 : DnaK protein

Lysates/proteins at 0.05 µg per lane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors