

Product datasheet

Recombinant Human hypothetical protein LOC254013
 ab166023

1 Image

Description

Product name	Recombinant Human hypothetical protein LOC254013	
Expression system	Wheat germ	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MALSLGWKAHRNHCGLLLQALRSSGLLLFPCGQCPW RGAGSFLDPEIKAF LEENTEVTSSGSLTPEIQLRLLTPRCKFWWERADLWP HSDPYWAIWPGG QALSRYLLDNPDVVRGKSVLDLGSGCGATAIAAKMSG ASRILANDIDPIA GMAITLNCELNRLNPFPIIQNILNLEQDKWDLVVLGDM FYDEDLADSLH QWLKCCFWTYRTRVLIGDPGRPQFSGHSIQHHLHKVV EYSLLESTRQENS GLTTSTVWGFQP	
Amino acids	1 to 262	
Tags	GST tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab166023** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00


Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Relevance

Hypothetical protein LOC254013 belongs to the UPF0516 family.

Images



ab166023 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human hypothetical protein LOC254013 (ab166023)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors