

Product datasheet

Recombinant Human IFI44L protein ab161220

1 Image

Description

Product name	Recombinant Human IFI44L protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>MVERCSRQGCTITMAYIDYNMIVAFMLGNYINLHESSTEPN DSLWFSLQK KNDTTEIETLLLNTAPKIIDEQLVCRLSKTDIFIICRDNKMLDK MITRN LKLRFYGHRQYLECEVFRVEGIKDNLDDIKRIIKAREHRNRL LADIRDYR PYADLVSEIRILLVGPVGSKGSSFFNSVKSIFHGHVGTQAV VGSDITSIT ERYRISVKDGNKSLPFMLCDTMGLDGAEGAGLCMDD IPHILKGCMPD RYQFNRPITPEHSTFITSPSLKDRIHCVAAYLDINSIDNLY SKMLAKV KQVHKEVLNCGIAYVALLTKVDDCSEVLQDNFLNMSRSM TSQSRVMNVHK MLGIPISNILMVGNYASDLELDPMKDILSALRQMLRAADD FLEDLPLE ETGAIERALQPCI</p>
Amino acids	1 to 413
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab161220** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot
	ELISA
Form	Liquid

Additional notes

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function

Exhibits a low antiviral activity against hepatitis C virus.


Sequence similarities

Belongs to the IFI44 family.

Cellular localization

Cytoplasm.

Images



ab161220 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human IFI44L protein
(ab161220)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors