

## Product datasheet

# Recombinant Human ITPKB protein ab127790

### Description

|                                   |                                 |
|-----------------------------------|---------------------------------|
| <b>Product name</b>               | Recombinant Human ITPKB protein |
| <b>Expression system</b>          | Escherichia coli                |
| <b>Accession</b>                  | <a href="#">P27987</a>          |
| <b>Protein length</b>             | Protein fragment                |
| <b>Animal free</b>                | No                              |
| <b>Nature</b>                     | Recombinant                     |
| <b>Species</b>                    | Human                           |
| <b>Predicted molecular weight</b> | 17 kDa                          |
| <b>Amino acids</b>                | 287 to 446                      |
| <b>Tags</b>                       | His-DHFR tag N-Terminus         |

### Specifications

Our [Abpromise guarantee](#) covers the use of **ab127790** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

|                     |             |
|---------------------|-------------|
| <b>Applications</b> | SDS-PAGE    |
| <b>Form</b>         | Lyophilised |

### Preparation and Storage

|                              |  |
|------------------------------|--|
| <b>Stability and Storage</b> | Shipped at 4°C. Store at -20°C.<br>Constituents: 0.32% Tris HCl, 0.58% Sodium chloride |
| <b>Reconstitution</b>        | Reconstitute in 100 µl of water.   |

### General Info

|                              |   |
|------------------------------|---|
| <b>Sequence similarities</b> | Belongs to the inositol phosphokinase (IPK) family. |
|------------------------------|---|

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors