

Product datasheet

Recombinant Human KPNA5 protein (Tagged) ab199581

Description

Product name	Recombinant Human KPNA5 protein (Tagged)
Purity	> 90 % SDS-PAGE. ab199581 was refolded using a unique temperature shift inclusion body refolding technology and chromatographically purified.
Expression system	Escherichia coli
Accession	O15131
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<pre> MASMTGGQQMGRGHHHHHHGNLYFQGGFEFDAMASP GKDNYRMKSYKNKAL NPQEMRRRRREEEGIQLRKQKREEQLFKRRNVYLPRND ESMLESPIQDPDI SSTVPIPEEEVVTDMVQMIFSNNADQQLTATQKFRKL LSKEPNPPIDQV IQKPGVVQRFVKFLERNENCTLQFEAAWALTNIASGTF LHTKVVIETGAV PIFIKLLNSEHEDVQEQA VWALGNIAGDNAECRDFVLN CEILPPLLELLT NSNRLTTTRNAVWALSNLCRGKNPPP NFSKVSPCLNV LSRLLFSSDPDVL ADVCWALSYS DGPNDKIQAVIDSGVCRRLLVELLMHN DYKVVSPALRAVG NIVTGDDIQTQVILNCSALPCLLHLLSSPKESIRKEACW TVSNITAGNRA QIQAVIDANIFPV LIEILQKA EFRTRKEAAWAITNATSGGT PEQIRYLVA LGCIKPLCDLLTVMDSKIVQVALNGLENILRLGEQESKQ NGIGINPYCAL IEEAYGLDKIEFLQSHENQEIQKAFDLIEHYFGVEEDD PSIVPQVDENQ QQFIFQQQEAPMDGFQL </pre>
Predicted molecular weight	64 kDa including tags

Amino acids 1 to 536

Additional sequence information Constructed with a small T7-His-TEV cleavage site Tag (29 amino acids) fusion at the N-terminal.

Specifications

Our [Abpromise guarantee](#) covers the use of **ab199581** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituent: 0.32% Tris HCl

Contains NaCl, KCl, EDTA, Sucrose, DTT.

General Info

Function Functions in nuclear protein import as an adapter protein for nuclear receptor KPNB1. Binds specifically and directly to substrates containing either a simple or bipartite NLS motif. Docking of the importin/substrate complex to the nuclear pore complex (NPC) is mediated by KPNB1 through binding to nucleoporin FxFG repeats and the complex is subsequently translocated through the pore by an energy requiring, Ran-dependent mechanism. At the nucleoplasmic side of the NPC, Ran binds to importin-beta and the three components separate and importin-alpha and -beta are re-exported from the nucleus to the cytoplasm where GTP hydrolysis releases Ran from importin. The directionality of nuclear import is thought to be conferred by an asymmetric distribution of the GTP- and GDP-bound forms of Ran between the cytoplasm and nucleus. Mediates nuclear import of STAT1 homodimers and STAT1/STAT2 heterodimers by recognizing non-classical NLSs of STAT1 and STAT2 through ARM repeats 8-9. Recognizes influenza A virus nucleoprotein through ARM repeat 7-9 In vitro, mediates the nuclear import of human cytomegalovirus UL84 by recognizing a non-classical NLS.

Tissue specificity Testis.

Sequence similarities Belongs to the importin alpha family.
Contains 10 ARM repeats.
Contains 1 IBB domain.

Domain Consists of an N-terminal hydrophilic region, a hydrophobic central region composed of 10 repeats, and a short hydrophilic C-terminus. The N-terminal hydrophilic region contains the importin beta binding domain (IBB domain), which is sufficient for binding importin beta and essential for nuclear protein import.
The IBB domain is thought to act as an intrasteric autoregulatory sequence by interacting with the internal autoinhibitory NLS. Binding of KPNB1 probably overlaps the internal NLS and contributes to a high affinity for cytoplasmic NLS-containing cargo substrates. After dissociation of the importin/substrate complex in the nucleus the internal autoinhibitory NLS contributes to a low affinity for nuclear NLS-containing proteins.

The major and minor NLS binding sites are mainly involved in recognition of simple or bipartite NLS motifs. Structurally located within in a helical surface groove they contain several conserved Trp and Asn residues of the corresponding third helices (H3) of ARM repeats which mainly contribute to binding.

Cellular localization

Cytoplasm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors