

Product datasheet

Recombinant Human L17 protein ab171590

1 Image

Description

<b>Product name</b>	Recombinant Human L17 protein	
<b>Purity</b>	> 85 % SDS-PAGE. ab171590 was purified using conventional chromatography.	
<b>Expression system</b>	Escherichia coli	
<b>Accession</b>	<a href="#">Q68CL5</a>	
<b>Protein length</b>	Full length protein	
<b>Animal free</b>	No	
<b>Nature</b>	Recombinant	
<b>Species</b>	Human	
<b>Sequence</b>	MGSSHHHHHHSSGLVPRGSHMGSMEEEEASSPGLGCSK PHLEKLTGITRI LESSPGVTEVTIIEKPPAERHMISSEWEQKNNCVMPELVKN FYLMTNGFHM TWSVKLDEHIPLGSMAINSISKLTQLTQSSMYSLPNAPTLA DLEDDTHE ASDDQPEKPHFDSRSVIFELDSCNGSGKVCLVYKSGKPA LAEDTEWFLD RALYWHFLDTFTAYRLLITHLGLPQWQYAFTSYGISPQAK QWFSMYKP ITYNTNLLTEETDSFVNKLDPSKVFKSKNKVIPKKKGPVQP AGGQKGPS GPSGPSTSSTKSSSGSGNPTRK	
<b>Predicted molecular weight</b>	36 kDa including tags	
<b>Amino acids</b>	1 to 300	
<b>Tags</b>	His tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab171590** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	SDS-PAGE Mass Spectrometry
<b>Mass spectrometry</b>	MALDI-TOF

<b>Form</b>	Liquid
<b>Additional notes</b>	This product was previously labelled as C18orf10.

## Preparation and Storage

---

<b>Stability and Storage</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 8.00 Constituents: 0.03% DTT, 0.32% Tris HCl, 0.06% EDTA, 50% Glycerol (glycerin, glycerine), 1.17% Sodium chloride
------------------------------	--

## General Info

---

<b>Cellular localization</b>	Cytoplasm > cytoskeleton.
------------------------------	---------------------------

---

## Images

---



15% SDS-PAGE analysis of ab171590 at 3µg.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors