

Product datasheet

Recombinant Human LCE3A protein (His tag) ab238231

1 Image

Description

Product name	Recombinant Human LCE3A protein (His tag)	
Purity	> 85 % SDS-PAGE.	
Expression system	Escherichia coli	
Accession	Q5TA76	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MSCQQNQQQCQPPPKCPAKSPAQCLPPASSSCAPSSG GCGPSSERSCCLS HHRCRRSHRCRCQSSNSCDRGSQQGGSSSCGHSSAG CC	
Predicted molecular weight	15 kDa including tags	
Amino acids	1 to 89	
Tags	His tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab238231** in the following tested applications.
 The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage


Stability and Storage	Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)
------------------------------	--

General Info

Relevance

LCE3A belongs to the LCE family. The function remains unknown.

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel analysis of ab238231.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors