

Product datasheet

Recombinant Human LRRC55 protein ab165918

1 Image

Overview

Product name	Recombinant Human LRRC55 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human
Sequence	MLRSPTFTDAGPRCSCLPVSQTLDSMDTVLMGSLQH CCCLLPKMGDTWAQ LPWPGPPHPAMLLISLLLAAGLMHSDAGTSCPVLCTC RNQVDCSSQRLF SVPPDLPM DTRNLSLAHNRTAVPPGYLTCYMELQVLD LHNNSLMELPRG LFLHAKRLAHL DLSYNNFSHV PADMFQEAHGLVHIDLS HNPWLRRVHPQA FQGLMQLRDL DLSYGGLAFLSLEALEGLPGLVTLQIGG NPWVCGCTMEPL LKWLRNRIQRCTADSQLAECRGPPEVEGAPLFSLTEE SFKACHLTLTDD YLFIAFVGFVVSIA SVATN FLLGITANCCHRWSKASEEE EI
Amino acids	1 to 341
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab165918** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot

Form	Liquid
-------------	--------

Additional notes

Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Sequence similarities

Contains 5 LRR (leucine-rich) repeats.


Contains 1 LRRCT domain.

Contains 1 LRRNT domain.

Cellular localization

Membrane.

Images



ab165918 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human LRRC55 protein

(ab165918)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors