

Product datasheet

Recombinant Human LZTR2 protein ab164789

1 Image

Description

Product name	Recombinant Human LZTR2 protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human

Sequence	<p>MELWAPQRLPQTRGKATAPSKDPDRGFRRDGHHRP VPHSWHNGERFHQWQ DNRGSPQPQQEPRADHQQQPHYASRPGDWHQPVSG VDYEGGYRNQLYSR PGYENYQSYQSPTMREEYAYGSYYYHGHPQWLQEER VPRQRSPYWHED YREQKYLDEHHYENQHSPFGTNSETHFQSNRNPKD SPASNSGQEWPGE LFPGSLLEAQKNKPSLASESNLLQQRESGLSSSSYE LSQYIRDAPERDD PPASAAWSPVQAEDVSSAGPKAPMKFYIPHPVVSFG PGGQLVRVGPSSPT DGQAALVELHSMEVILNDSEEQEEMRSFSGPLIREDV HKVDIMTFCQQA AQCKSETLGSRDSALLWQLLVLLCRQNGSMVGS DIA ELLMQDCKKLEKY KRQPPVANLINLTDEDWPVLSSGTPNLLTGEIPPSVET PAQVEKFTRLL YYGRKKEALEWAMKNHLWGHALFLSSKMDPQTYSWV MSGFTSTLALNDPL QTLFQLMSGRIQAATCCGEKQWGDWRPHLAVILSNQ AGDPELYQRAVA IGDTLAGKGLVEAAHFCYMAHVDPFGHYTVKTDHLVLL GSSHRYATWEKG NSKDIFQGTVLALVGFYGSSFHFLM</p>
-----------------	---

Amino acids	1 to 625
Tags	GST tag N-Terminus

Specifications

Specifications

Our [Abpromise guarantee](#) covers the use of **ab164789** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot
Form	Liquid
Additional notes	


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Function	Required for secretory cargo traffic from the endoplasmic reticulum to the Golgi apparatus and for normal transitional endoplasmic reticulum (tER) organization.
Tissue specificity	Ubiquitous.
Sequence similarities	Belongs to the SEC16 family.
Cellular localization	Endoplasmic reticulum membrane. Golgi apparatus membrane.

Images



ab164789 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors