

Product datasheet

Recombinant Human MAP3K15/ASK3 protein ab166430

1 Image

Overview

Product name	Recombinant Human MAP3K15/ASK3 protein
Protein length	Protein fragment

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human
Sequence	LYHLQLKLKSNCTENPAGPYGQRTDKELIDWLRQLGA DAKTIEKIVEEG YTLSDILNEITKEDLRYLRRLRGGLLCRLWSAVSQYRRAQ EASETKD
Amino acids	691 to 786
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab166430** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml. This product was previously labelled as MAP3K15

Preparation and Storage


Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------------------------	---

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function	May function in a signal transduction pathway that is activated by various cell stresses and leads to apoptosis.
Tissue specificity	Isoform 2 and isoform 3 are widely expressed. Isoform 2 highest levels are observed in fetal brain, and isoform 3 highest levels in pancreas, peripheral blood leukocytes, fetal brain and spleen.
Sequence similarities	Belongs to the protein kinase superfamily. STE Ser/Thr protein kinase family. MAP kinase kinase kinase subfamily. Contains 1 protein kinase domain.

Images



ab166430 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human MAP3K15/ASK3 protein (ab166430)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors