

Product datasheet

Recombinant Human MOS protein ab125632

1 Image

Description

| | |
|-----------------------------------|---------------------------------|
| Product name | Recombinant Human MOS protein |
| Expression system | Baculovirus infected Sf21 cells |
| Accession | P00540 |
| Protein length | Full length protein |
| Animal free | No |
| Nature | Recombinant |
| Species | Human |
| Predicted molecular weight | 66 kDa including tags |
| Amino acids | 1 to 346 |

Specifications

Our [Abpromise guarantee](#) covers the use of **ab125632** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|---------------------|--------------------------|
| Applications | Western blot SDS-PAGE |
| Form | Liquid |

Preparation and Storage

| | |
|------------------------------|---|
| Stability and Storage | Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.31% Glutathione, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA, 25% Glycerol, 0.88% Sodium chloride |
|------------------------------|---|

General Info

| | |
|------------------------------|--|
| Tissue specificity | Expressed specifically in testis during spermatogenesis. |
| Sequence similarities | Belongs to the protein kinase superfamily. Ser/Thr protein kinase family. Contains 1 protein kinase domain. |

Images



SDS-PAGE of ab125632

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors