

Product datasheet

Recombinant Human NeuroD2 protein ab152568

1 Image

Description

<b>Product name</b>	Recombinant Human NeuroD2 protein
<b>Expression system</b>	Wheat germ
<b>Accession</b>	<a href="#">Q15784</a>
<b>Protein length</b>	Full length protein
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Sequence</b>	<pre> MLTRLFSEPGLLSDVPKFASWGDGEDDEPRSDKGDA PPPPPPAPGPGAPG PARAAKPVPLRGEEGTEATLAEVKEEGELGGEEEE EEEEGLDEAEGER PKKGGPKKRKMTKARLERSKLRQKANARERNRMD LNAALDNLRKVVPC YSKTQKLSKIETLRLAKNYWALSEILRSGKRPDLVSYV QTLCKGLSQPT TNLVAGCLQLNSRNFLTEQGADGAGRFHGSGGPFAM HPYPYPCSRLAGAQ CQAAGGLGGGAAHALRTHGYCAAYETLYAAAGGGGA SPDYNSEYEGPLS PPLCLNGNFSKQDSSPDHEKSYHYSMHYSALPGSRP TGHGLVFGSSAVR GGVHSENLISYDMHLHHDGRGPMYEELNAFFHN                     </pre>
<b>Predicted molecular weight</b>	68 kDa including tags
<b>Amino acids</b>	1 to 382

Specifications

Our [Abpromise guarantee](#) covers the use of **ab152568** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	Western blot
	ELISA
	SDS-PAGE

<b>Form</b>	Liquid
<b>Additional notes</b>	Protein concentration is above or equal to 0.05 µg/µl.

## Preparation and Storage

---

<b>Stability and Storage</b>	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

## General Info


---

<b>Function</b>	Appears to mediate neuronal differentiation.
<b>Sequence similarities</b>	Contains 1 basic helix-loop-helix (bHLH) domain.
<b>Cellular localization</b>	Nucleus.

---

## Images

---



12.5% SDS-PAGE analysis of ab152568 stained with Coomassie Blue.

SDS-PAGE - Recombinant Human NeuroD2 protein (ab152568)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors