

Product datasheet

Recombinant human NGF protein (Animal Free) ab179616

3 References

Description

Product name	Recombinant human NGF protein (Animal Free)	
Biological activity	The activity is determined by the ability to stimulate proliferation of TF-1 cells and is typically less than 1 ng/mL.	
Purity	> 97 % SDS-PAGE. assessed by HPLC, Reducing and Non-reducing SDS-PAGE, UV spectroscopy at 280 nm. Note: This product is produced with no animal-derived raw products, animal free equipment and animal free protocols.	
Endotoxin level	<=1.000 Eu/µg	
Expression system	Escherichia coli	
Accession	P01138	
Protein length	Protein fragment	
Animal free	Yes	
Nature	Recombinant	
Species	Human	
Sequence	SSSHPIFHRGEFSVCDSVSVWVGDKTTATDIKGKEVMVL GEVNINNSVFK QYFFETKCRD PNPVDSGCRGIDSKHWNSYCTTTHTFVKALTM DGKQAA WRFIRIDTACVCVLSRKAVRRA	
Predicted molecular weight	14 kDa	
Amino acids	122 to 241	
Additional sequence information	Mature form.	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab179616** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	HPLC
	SDS-PAGE
	Functional Studies

Form Lyophilized

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at -20°C. For long term storage it is recommended to add a carrier protein on reconstitution (0.1% HSA or BSA).

This product is an active protein and may elicit a biological response in vivo, handle with caution.

Reconstitution Centrifuge vial before opening. When reconstituting the product, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution. It is recommended to reconstitute the lyophilized product with sterile water at a concentration of 0.1 mg/mL, which can be further diluted into other aqueous solutions. Reconstituted material should be aliquoted and frozen at -20°C

General Info

Function Nerve growth factor is important for the development and maintenance of the sympathetic and sensory nervous systems. Extracellular ligand for the NTRK1 and NGFR receptors, activates cellular signaling cascades through those receptor tyrosine kinase to regulate neuronal proliferation, differentiation and survival. Inhibits metalloproteinase dependent proteolysis of platelet glycoprotein VI (PubMed:20164177).

Involvement in disease Neuropathy, hereditary sensory and autonomic, 5

Sequence similarities Belongs to the NGF-beta family.

Cellular localization Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors