

Product datasheet

Recombinant Human NKG2C protein (denatured)
ab131674

1 Image

Description

Product name	Recombinant Human NKG2C protein (denatured)
Purity	> 85 % SDS-PAGE.
Expression system	Escherichia coli
Accession	P26717
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMIPFLEQ NNFSPNTRTQ KARHCGHCPE EWITYNSCY YIGKERRTWE ELLACTSKN SLLSIDNEE EMKFLASILP SSWIGVFRNS SHHPWWTING LAFKHKIKDS DNAELNCAVL QVNRLKSAQC GSSMIYHCKH KL
Predicted molecular weight	18 kDa including tags
Amino acids	94 to 231
Tags	His tag N-Terminus
Description	Recombinant Human NKG2C protein

Specifications

Our [Abpromise guarantee](#) covers the use of **ab131674** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 2.4% Urea, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine)

General Info

Function	Plays a role as a receptor for the recognition of MHC class I HLA-E molecules by NK cells and some cytotoxic T-cells.
Tissue specificity	Natural killer cells.
Sequence similarities	Contains 1 C-type lectin domain.
Cellular localization	Membrane.

Images



15% SDS-PAGE analysis of 3 µg ab131674.

SDS-PAGE - Recombinant Human NKG2C protein
(denatured) (ab131674)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors