

Product datasheet

Recombinant Human PDCL2 protein ab127854

Description

Product name	Recombinant Human PDCL2 protein
Purity	= 90 % Purified via His tag.
Expression system	Escherichia coli
Accession	Q8N4E4
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Amino acids	5 to 237
Tags	His-DHFR tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab127854** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Lyophilized

Preparation and Storage

Stability and Storage	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
Reconstitution	Reconstitute with water to desired concentration.

General Info

Function	May play a role in germ cell maturation.
Sequence similarities	Belongs to the phosphducin family.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors