

Product datasheet

Recombinant Human PIK3IP1 protein ab185388

Description

Product name	Recombinant Human PIK3IP1 protein	
Purity	> 95 % SDS-PAGE. Determined by SEC-HPLC and reducing SDS-PAGE.	
Endotoxin level	< 1.000 Eu/μg	
Expression system	Mammalian	
Accession	Q96FE7	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	SGGCFWDNGHLYREDQTSPAPGLRCLNWLDAQSGL ASAPVSGAGNHSYCR NPDEDPRGPWCYVSGEAGVPEKRPCEDLRCPETTS QALPAFTTEIQEASE GPGADEVQVFAPANALPARSEAAAVQPVIGISQRVRM NSKEKKDLGTVDH HHHHH	
Predicted molecular weight	17 kDa including tags	
Amino acids	22 to 168	
Tags	His tag C-Terminus	
Additional sequence information	Extracellular domain.	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab185388** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE HPLC
Form	Lyophilised

Preparation and Storage

Stability and Storage Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term. Please see notes section.
pH: 8.00
Constituents: 0.79% Tris HCl, 0.31% Glutathione

Reconstitution Dissolve the lyophilized protein in 3X PBS. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Reconstituted protein should be aliquoted and stored at < -20 degrees.

General Info

Sequence similarities Contains 1 krigle domain.

Cellular localization Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors