

Product datasheet

Recombinant Human PPP3CB protein ab159176

1 Image

Overview

Product name Recombinant Human PPP3CB protein
Protein length Protein fragment

Description

Nature Recombinant
Source Wheat germ
Amino Acid Sequence
Species Human
Sequence LTPTGMLPSGVLGGRRQLQSATVEAIEAEKAIRGFSP
 PHRICSFEEAKG
 LDRINERMPPRKDAVQQDGFNSLNTAHATENHGTGNH
 TAQ
Amino acids 435 to 524
Tags proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab159176** in the following tested applications.
 The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Western blot
 ELISA
Form Liquid
Additional notes Protein concentration is above or equal to 0.05 mg/ml.


Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
 pH: 8.00
 Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function	Calcium-dependent, calmodulin-stimulated protein phosphatase. This subunit may have a role in the calmodulin activation of calcineurin.
Sequence similarities	Belongs to the PPP phosphatase family. PP-2B subfamily.

Images



ab159176 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human PPP3CB protein
(ab159176)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors