

Product datasheet

Recombinant Human PSF3 protein ab163880

1 Image

Overview

Product name	Recombinant Human PSF3 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human
Sequence	MSEAYFRVESGALGPEENFLSLDDILMSHEKLPVRTET AMPRLGAFFLER SAGAETDNAVPQGSKLELPLWLAKGLFDNKRRILSVE LPKIYQEGWRTVF SADPNVVDLHKMGPHFYGFSGQLLHFDSPENADISQS LLQTFIGRFRRIM DSSQNAYSEDTSALVARLDEMERGLFQTGQKGLNDF QCWEKGQASQITAS NLVQNYKKRKFTDMED
Amino acids	1 to 216
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab163880** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------------------------	---


pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function	The GINS complex plays an essential role in the initiation of DNA replication, and progression of DNA replication forks. GINS complex seems to bind preferentially to single-stranded DNA.
Sequence similarities	Belongs to the GINS3/PSF3 family.
Cellular localization	Nucleus.

Images



ab163880 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human PSF3 protein
(ab163880)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors