

Product datasheet

Recombinant Human RAB31 protein ab124565

1 Image

Description

Product name	Recombinant Human RAB31 protein	
Purity	> 90 % SDS-PAGE. ab124565 was purified by conventional chromatography techniques.	
Expression system	Escherichia coli	
Accession	Q13636	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHLHKFLWDTA GQERFHSLAP MMYRGSAAAV IYDITKQDS FYTLKKWVKE LKEHGPENIV MAIAGNKCDL SDIREVPLKD AKEYAESIGA IVVETSAKNA INIEELFQGI SRQIPPLDPH ENGNGTIKV EKPTMQASRR C	
Predicted molecular weight	26 kDa including tags	
Amino acids	1 to 195	
Tags	His tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab124565** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------------------------	--


pH: 8.00

Constituents: 0.32% Tris HCl, 30% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

General Info

Tissue specificity	Highest expression in placenta and brain with lower levels in heart and lung. Not detected in liver, skeletal muscle, kidney or pancreas.
Sequence similarities	Belongs to the small GTPase superfamily. Rab family.
Cellular localization	Cell membrane.

Images



ab124565 (3µg) visualised by SDS-PAGE (15%).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors