

Product datasheet

Recombinant Human RAB8B protein (denatured)
ab167867

1 Image

Overview

Product name	Recombinant Human RAB8B protein (denatured)
Protein length	Full length protein
Description	Recombinant Human RAB8B protein

Description

Nature	Recombinant
Source	Escherichia coli
Amino Acid Sequence	
Accession	Q92930
Species	Human
Sequence	<p>MGSSHHHHHH SSGLVPRGSH MGSMAKTYDY LFKLLIGDS GVGKTCLLFR FSEDAFNNTF ISTIGIDFKI RTIELDGKKI KLQWDTAGQ ERFRTITTAY YRGAMGIMLV YDITNEKSFN NIKNWIRNIE EHASSDVERM ILGNKCDMND KRQVSKERGE KLAIDYGIKF LETSAKSSAN VEEAFFTLAR DIMTKLNRKM NDSNSAGAGG PVKITENRSK KTSFFRC</p>
Molecular weight	26 kDa including tags
Amino acids	1 to 204
Tags	His tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab167867** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Purity	>90% by SDS-PAGE.

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00
Constituents: 0.32% Tris-HCl buffer, 2.4% Urea, 10% Glycerol


General Info

Function May be involved in vesicular trafficking and neurotransmitter release.

Sequence similarities Belongs to the small GTPase superfamily. Rab family.

Cellular localization Cell membrane.

Images



15% SDS-PAGE analysis of 3µg ab167867.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors