

Product datasheet

Recombinant Human RPL7A protein (denatured)
 ab187474

1 Image

Description

| | |
|-----------------------------------|--|
| Product name | Recombinant Human RPL7A protein (denatured) |
| Purity | > 85 % SDS-PAGE. |
| Expression system | Escherichia coli |
| Accession | P62424 |
| Protein length | Full length protein |
| Animal free | No |
| Nature | Recombinant |
| Species | Human |
| Sequence | <pre> MGSSHHHHHH SSGLVPRGSH MGSMPKGKKA KGKKVAPAPA VVKKQEAKKV VNPLFEKRPK NFGIGQDIQP KRDLTRFVKW PRYIRLQRQR AILYKRLKVP PAINQFTQAL DRQTATQLLK LAHKYRPETK QEKKQRLAR AEKKAAGKGD VPTKRPPVLR AGVNTVTTLV ENKKAQLVVI AHDVDPIELV VFLPALCRKM GVPYCIKGG ARLGRLVHRK TCTTVAFTQV NSEDKGALAK LVEAIRTNYN DRYDEIRRHW GGNVLGPKSV ARIAKLEKAK AKELATKLG </pre> |
| Predicted molecular weight | 32 kDa including tags |
| Amino acids | 1 to 266 |
| Tags | His tag N-Terminus |
| Description | Recombinant Human RPL7A protein |

Specifications

Our [Abpromise guarantee](#) covers the use of **ab187474** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8

Constituents: 0.32% Tris HCl, 10% Glycerol

General Info


Involvement in disease

Note=Chromosomal recombination involving RPL7A activates the receptor kinase domain of the TRK oncogene.

Sequence similarities

Belongs to the ribosomal protein L7Ae family.

Images



ab187474 (3 µg) on 15 % SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors