

Product datasheet

Recombinant Human RPS2 protein (denatured)
 ab171474

1 Image

Description

Product name	Recombinant Human RPS2 protein (denatured)
Purity	> 85 % SDS-PAGE.
Expression system	Escherichia coli
Accession	P15880
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<pre> MGSSHHHHHH SSGLVPRGSH MGSMADDAGA AGGPGGPGGP GMGNRGGFRG GFGSGIRGRG RGRGRGRGRG RGARGGKAED KEWMPVTKLG RLVKDMKIKS LEEMLFSLP IKESEIDFF LGASLKDEVL KIMPVQKQTR AGQRTRFKAF VAIGDYNGHV GLGVKCSKEV ATAIRGAIL AKLSVPVRR GYWGNKIGKP HTVPCKVTGR CGSVLVRLIP APRGTGMSA PVPKLLMMA GIDDCYTSAR GCTATLGNFA KATFDAISK YSYLTPDLWK ETVFTKSPYQ EFTDHLVKTH TRVSVQRTQA PAVATT </pre>
Predicted molecular weight	34 kDa including tags
Amino acids	1 to 293
Tags	His tag N-Terminus
Description	Recombinant Human RPS2 protein

Specifications

Our [Abpromise guarantee](#) covers the use of **ab171474** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 2.4% Urea, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine)


General Info

Sequence similarities

Belongs to the ribosomal protein S5P family.

Contains 1 S5 DRBM domain.

Images



15% SDS-PAGE analysis of 3µg ab171474.

SDS-PAGE - Recombinant Human RPS2 protein
(denatured) (ab171474)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors