

Product datasheet

Recombinant Human RYK protein ab139623

1 Image

Description

Product name	Recombinant Human RYK protein
Purity	> 90 % SDS-PAGE.
Expression system	Baculovirus infected Sf9 cells
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>MKRIELDDISISASSSSQGLSQPSTQTTQYL RADTPN NATPI TSSLGYPTL RIEKNDLR SVTLLEAKGKVKDIAISRERITLKDVLQEGTFGR I FHGILID EKDPNKEKQAFVKTVDQASEIQVTMMLTESCKLRGLHH RNLLPITHVCI EEGEKPMVILPYMNWGNLKLFLRQCKLVEANNPQAISQQ DLVHMAIQIAC GMSYLARREVIHKDLAARNCVIDDTLQVKITDNALSRDLFP MDYHCLGDN ENRPVRWMALESLVNNEFSSASDVWAFGVTLWELMTLG QTPYVDIDPFEM AAYLKDGYRIAQPINCPDELFAVMACCWALDPEERPKFQ QLVQCLTEFHA ALGAYV</p>
Predicted molecular weight	66 kDa including tags
Amino acids	255 to 610

Specifications

Our [Abpromise guarantee](#) covers the use of **ab139623** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 7.50

Constituents: 0.31% Glutathione, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA, 25% Glycerol (glycerin, glycerine), 0.29% Sodium chloride

General Info

Function

May be a coreceptor along with FZD8 of Wnt proteins, such as WNT1, WNT3, WNT3A and WNT5A. Involved in neuron differentiation, axon guidance, corpus callosum establishment and neurite outgrowth. In response to WNT3 stimulation, receptor C-terminal cleavage occurs in its transmembrane region and allows the C-terminal intracellular product to translocate from the cytoplasm to the nucleus where it plays a crucial role in neuronal development.

Tissue specificity

Observed in all the tissues examined.

Sequence similarities

Belongs to the protein kinase superfamily. Tyr protein kinase family.

Contains 1 protein kinase domain.

Contains 1 WIF domain.

Domain

The extracellular WIF domain is responsible for Wnt binding.


Post-translational modifications

Proteolytically cleaved, in part by presenilin, in response to WNT3 stimulation. Cleavage occurs during neuronal differentiation.

Cellular localization

Membrane. Nucleus. Cytoplasm. In cells that have undergone neuronal differentiation, the C-terminal cleaved part is translocated from the cytoplasm to the nucleus.

Images



SDS PAGE analysis of ab139623

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors