

Product datasheet

Recombinant Human SLFN11 protein ab164865

1 Image

Description

Product name	Recombinant Human SLFN11 protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human

Sequence	<p>MEANQCPLVVEPSYPDLVINVGCVTLGEENRKKLQKIQ RDQEKERVMAAA CALLNSGGGVIRMAKKVEHPVEMGLDLEQSLRELIQS SDLQAFFETKQQG RCFYIFVKSWSGGPFPEDRSVKPRCLSLSSSLYRRSE TSVRSMDSREAFCA FLKTKRKPKILEEGPFHKIHKGVYQELPNSDPADPNSD PADLIFQKDYLE YGEILPFESQLVEFKQFSTKHFQEYVKRTIPEYVPAFA NTGGGYLFIGV DDKSREVLGCAKENVDPDSLRRKIEQAMKLPVHFC QPQRPIFTLKV DVLKRGELYGYACMIRVNPFCVAVSEAPNSWVMDK YVCSLTTEKWVGM MTDTPDLLQLSEDFECQLSLSSGPPLSRPVYSKKGL EHKKELQQLFSV PPGYLRYTPESLWRDLISEHRGLEELINKMQPFFRGILI FSRSWAVDLN LQEKPGVICDALLIAQNSTPILYTLREQDAEGQDYCTR AFTLKQKLVN MGGYTGKVCVRAKVLCLSPESAEALEAAVSPMDYP ASYSLAGTQHMEAL LQSLVMLLGFRSLLSDQLGCEVLNLLTAQQYEIFSRSL RKNRELFVHGL PGSGKTIMAMKIMEKIRNVFHCEAHRILYVCENQPLRNF ISDRNICRAET RKTFLRENFEBHQHVIDEAQNFRTEDGDWYGGKAKSITR RAKGGPGILWI FLDYFQTSHLDCSGLPPLSDQYPREELTRVIRNADPIA</p>
-----------------	--

KYLQKEMQVIRS
NPSFNIPTGCLEVFPEAEWSQGVQGTLRRIKKYLTVEQI
MTCVADTCRRFF
DRGYSPKDVAVLVSTAKEVEHYKYELLKAMRKKRVVQ
LSDACDMLGDHV
LDSVRRFSGLERSVFGIHPRTADPAILPNVLI CLASRAK
QHL YIFPWGG H

Amino acids 1 to 901
Tags GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab164865** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Western blot
ELISA
Form Liquid
Additional notes Protein concentration is above or equal to 0.05 mg/ml.


Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Sequence similarities Belongs to the Schlafen family.

Images



ab164865 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors