

Product datasheet

Recombinant Human SNURF protein ab177617

1 Image

Description

Product name	Recombinant Human SNURF protein	
Purity	> 90 % SDS-PAGE. ab177617 was purified using conventional chromatography techniques.	
Expression system	Escherichia coli	
Accession	Q9Y675	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMERARDR LHLRRTTEQH VPEVEVQVKR RRTASLSNQE CQLYPRRSQQ QQVPVDFQA ELRQAFLAET PRGG	
Predicted molecular weight	11 kDa including tags	
Amino acids	1 to 71	
Tags	His tag N-Terminus	
Additional sequence information	NP_073715.1	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab177617** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.32% Tris HCl, 0.88% Sodium chloride, 40% Glycerol (glycerin, glycerine), 0.02% DTT

General Info

Tissue specificity	Expressed in heart, skeletal muscle and lymphoblasts (at protein level). Expressed in brain, pancreas, heart, liver, lung, kidney and skeletal muscle.
Sequence similarities	Belongs to the SNURF family.
Cellular localization	Nucleus.

Images



15% SDS-PAGE analysis of ab177617 (3µg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors