

Product datasheet

Recombinant Human SPINK4 protein ab152050

Description

Product name	Recombinant Human SPINK4 protein	
Purity	> 95 % SDS-PAGE. Purity is greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE. ab152050 is 0.2 µM filtered.	
Endotoxin level	< 1.000 Eu/µg	
Expression system	HEK 293 cells	
Accession	O60575	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	GKLPFSRMPICEHMHVESPTCSQMSNLVCGTDGLTYTNEC QLCLARIKTKQ DIQIMKDGKCVDHSHHHH	
Predicted molecular weight	8 kDa including tags	
Amino acids	27 to 86	
Tags	His tag C-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab152050** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE HPLC
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.4 Constituents: 0.39% MES, 0.05% Brij, 0.02% Calcium chloride, 0.02% DTT, 10% Glycerol (glycerin, glycerine), 0.88% Sodium chloride
------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

General Info

Sequence similarities	Contains 1 Kazal-like domain.
Cellular localization	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors