

Product datasheet

Recombinant Human STAC3 protein ab166011

1 Image

Description

Product name	Recombinant Human STAC3 protein	
Expression system	Wheat germ	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	<p>MELPPEPQANGEAVGAGGGPIYYIEEEEEEEEEEEEE</p> <p>PPPEPPKLVNDKP</p> <p>HKFKDHFFKKPKFCDVCARMMLNKFGLRCKNCKTNI</p> <p>HEHCQSYVEMQR</p> <p>CFGKIPPGFHraySSPLYSNQQYACVKDLsAANRNDP</p> <p>VFETLRTGVIMAN</p> <p>KERKKGQADKKNPVAAMMEEEPESARPEEGKPQDG</p> <p>NPEGDKKAekKTPDD</p> <p>KHKQPGFQQSHYFVALYRFKALEKDDLDfPPGEKITVI</p> <p>DDSNEEWWRGKI</p> <p>GEKVGFFPPNFIRVRAGERVHRVTRSFVGNREIGQITL</p> <p>KKDQMVQKGD EAGGYVKVYtGRKVGLFPTDFLEEI</p>	
Amino acids	1 to 325	
Tags	GST tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab166011** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl


General Info

Sequence similarities

Contains 1 phorbol-ester/DAG-type zinc finger.

Contains 2 SH3 domains.

Images



ab166011 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors