

Product datasheet

Recombinant Human TAF1L protein ab196115

1 Image

Description

Product name	Recombinant Human TAF1L protein	
Purity	>= 82 % SDS-PAGE.	
Expression system	Escherichia coli	
Accession	Q8IZX4	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	TDPMVTLSSILESIIINDMRDLPNTHPFHTPVNAKVVKDYKII TRPMDLQ TLRENVKCLYPSREEFREHLELMKNSATYNGPKHSLTQI SQSMLDLCD EKLKEKEDKLARLEKAINP	
Predicted molecular weight	41 kDa including tags	
Amino acids	1398 to 1516	
Tags	GST tag N-Terminus	
Additional sequence information	Bromodomain 1 (NM_153809).	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab196115** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.79% Tris HCl, 0.8% Sodium chloride, 0.02% Potassium chloride, 0.61% Glutathione, 0.09% EDTA, 10% Glycerol, 0.03% DTT
------------------------------	--

General Info

Function	May act as a functional substitute for TAF1/TAFII250 during male meiosis, when sex chromosomes are transcriptionally silenced.
Tissue specificity	Testis specific, expressed apparently in germ cells.
Sequence similarities	Belongs to the TAF1 family. Contains 2 bromo domains.
Cellular localization	Nucleus.

Images



4-20% SDS-PAGE analysis of 1 μ g ab196115 with Coomassie staining.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors