

Product datasheet

Recombinant Human TBRG4 protein ab132734

1 Image

Description

Product name	Recombinant Human TBRG4 protein
Expression system	Wheat germ
Accession	Q969Z0
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human

Sequence

```

MAAHLVKRCTCLLREARQAPAMAPVGRLRLAWVAHKTL
TSSATSPISHL
PGSLMEPVEKERASTPYIEKQVDHLIKKATRPEELLELLGG
SHDLDSNQA
AMVLIRLSHLLSEKPEDKGLLIQDAHFHQLLCLLNSQIASV
WHGTLKLL
GSLYALGIPKASKELQSVEQEVRRMRKLYKHLAFLAES
CATLSQEQHS
QELLAELLTHLERRWTEIEDSHTLVTVMMKVGHLSPLMN
RLEDKCLELV
EHFGPNELRKVLVMLAAQSRRSVPLLRRAISYHLVQKPFSL
TKDVLLDVAY
AYGKLSFHQTQVSQRLATDLLSLMPSLTSGEVAHCAKSF
ALLKWLSLPLF
EAFAQHVLNRAQDITLPHLCSVLLAFARLNHFDPQEDQFF
SLVHEKLGSE
LPGLEPALQVDLVWALCVLQQAREAEQAVLHPEFHIQFL
GGKSQKDQNT
FQKLLHINATALLEYPEYSGPLLPAVAVPGPSALDRKVTP
LQKELQETL
KGLLGSADKGSLEVATQYGWVLDAEVLLDSDGEFLPVRD
FVAPHLAQPTG
SQSPPPGSKRLAFLRWEFPNFNSRSKDLLGRFVLARRHI
VAAGFLMDVP
FYEWLELKSEWQKGAYLKDKMRKAVAEELAK

```

Predicted molecular weight	97 kDa including tags
Amino acids	1 to 631

Specifications

Our [Abpromise guarantee](#) covers the use of **ab132734** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications ELISA
SDS-PAGE
Western blot

Form Liquid

Additional notes

Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.31% Glutathione, 0.79% Tris HCl


General Info

Function May play a role in cell cycle progression.

Tissue specificity Ubiquitously expressed.

Sequence similarities Belongs to the FAST kinase family.
Contains 1 RAP domain.

Images



12.5% SDS-PAGE analysis of ab132734 stained with Coomassie Blue.

SDS-PAGE - Recombinant Human TBG4 protein
(ab132734)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors