

Product datasheet

Recombinant Human Tissue type Plasminogen Activator protein ab185834

Description

Product name	Recombinant Human Tissue type Plasminogen Activator protein	
Purity	> 95 % SDS-PAGE. Purity is determined by SEC-HPLC and reducing SDS-PAGE.	
Endotoxin level	< 1.000 Eu/μg	
Expression system	Mammalian	
Accession	P00750	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	SYQVICRDEKTQMIYQQHQSWLRPVLRSNRVEYCWCN SGRAQCHSVPVKS CSEPRCFNNGGTCQQALYFSDVCQCPEGFAGKCCEI DTRATCYEDQGISY RGTWSTAESGAECTNWNSSALAQKPYSGRRPDAIRL GLGNHNYCRNPDRD SKPWCYVFKAGKYSSEFCSTPACSEGNSDCYFGNGS AYRGTHSLTESGAS CLPWNSMILIGKVYTAQNPQAALGLGKHNYCRNPDG DAKPWCHVLKNRR LTWEYCDVPSCSTCGLRQYSQPQFRIKGGLFADIASH PWQAAIFAKHRRS PGERFLCGGILISSCWILSAAHCFQERFPPHHLTIVLGR TYRVVPGEEEEQ KFEVEKYVHKEFDDDDTYDNDIALQLKSDSSRCAQES SVVRTVCLPPAD LQLPDWTECELSGYGKHEALSPFYSERLKEAHVRLYP SSRCTSQHLLNRT VTDNMLCAGDTRSGGPQANLHDACQGDSGGPLVCL NDGRMTLVGIISWGL GCGQKDVPGVYTKVTNYLDWIRDNMRPVDHHHHHH	
Predicted molecular weight	61 kDa including tags	
Amino acids	36 to 562	

Tags	His tag C-Terminus
Additional sequence information	Full length mature protein

Specifications

Our [Abpromise guarantee](#) covers the use of **ab185834** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
	HPLC
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 5.50 Constituents: 30% Glycerol, 0.88% Sodium chloride, 0.392% MES, 0.0022% Calcium chloride
------------------------------	---

General Info

Function	Converts the abundant, but inactive, zymogen plasminogen to plasmin by hydrolyzing a single Arg-Val bond in plasminogen. By controlling plasmin-mediated proteolysis, it plays an important role in tissue remodeling and degradation, in cell migration and many other physiopathological events. Play a direct role in facilitating neuronal migration.
Tissue specificity	Synthesized in numerous tissues (including tumors) and secreted into most extracellular body fluids, such as plasma, uterine fluid, saliva, gingival crevicular fluid, tears, seminal fluid, and milk.
Involvement in disease	Note=Increased activity of TPA results in increased fibrinolysis of fibrin blood clots that is associated with excessive bleeding. Defective release of TPA results in hypofibrinolysis that can lead to thrombosis or embolism.
Sequence similarities	Belongs to the peptidase S1 family. Contains 1 EGF-like domain. Contains 1 fibronectin type-I domain. Contains 2 kringle domains. Contains 1 peptidase S1 domain.
Domain	Both FN1 and one of the kringle domains are required for binding to fibrin. Both FN1 and EGF-like domains are important for binding to LRP1. The FN1 domain mediates binding to annexin A2. The second kringle domain is implicated in binding to cytokeratin-8 and to the endothelial cell surface binding site.
Post-translational modifications	The single chain, almost fully active enzyme, can be further processed into a two-chain fully active form by a cleavage after Arg-310 catalyzed by plasmin, tissue kallikrein or factor Xa. Differential cell-specific N-linked glycosylation gives rise to two glycoforms, type I (glycosylated at Asn-219) and type II (not glycosylated at Asn-219). The single chain type I glycoform is less readily converted into the two-chain form by plasmin, and the two-chain type I glycoform has a lower activity than the two-chain type II glycoform in the presence of fibrin. N-glycosylation of Asn-152; the bound oligomannosidic glycan is involved in the interaction with

the mannose receptor.

Characterization of O-linked glycan was studied in Bowes melanoma cell line.

Cellular localization

Secreted > extracellular space.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors