

Product datasheet

Recombinant Human TRPC1 protein ab152768

1 Image

Description

Product name	Recombinant Human TRPC1 protein	
Expression system	Wheat germ	
Accession	P48995-2	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	AHNKFHDFADRKDWDAFHPTLVAEGLFAFANVLSYLRLF FMYTTSSILGP LQISMGQMLQDFGK	
Predicted molecular weight	33 kDa including tags	
Amino acids	442 to 505	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab152768** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot
	SDS-PAGE
	ELISA
Form	Liquid

Additional notes


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------

General Info

Function	Thought to form a receptor-activated non-selective calcium permeant cation channel. Probably is operated by a phosphatidylinositol second messenger system activated by receptor tyrosine kinases or G-protein coupled receptors. Seems to be also activated by intracellular calcium store depletion.
Tissue specificity	Seems to be ubiquitous.
Sequence similarities	Belongs to the transient receptor (TC 1.A.4) family. STrpC subfamily. TRPC1 sub-subfamily. Contains 4 ANK repeats.
Post-translational modifications	Activation of PRKCA induces phosphorylation of TRPC1 and subsequent Ca ²⁺ entry into cells.
Cellular localization	Membrane.

Images



12.5% SDS-PAGE analysis of ab152768 stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors