

Product datasheet

Recombinant Human Tspan-8 protein ab159720

1 Image

Description

| | |
|--------------------------|---|
| Product name | Recombinant Human Tspan-8 protein |
| Expression system | Wheat germ |
| Protein length | Full length protein |
| Animal free | No |
| Nature | Recombinant |
| Species | Human |
| Sequence | MAGVSACIKYSMFTFNFLFWLCGILILALAWVRISNDSQ AIFGSEDVGS SSYVAVDILIAVGAIIMLGFLACCGAIKESRCMLLLFFIGL LLILLQV ATGILGAVFKSKSDRVNETLYENTKLLSATGESEKQFQ EAIMVQEEFK CCGLVNGAADWGNNFQHYPELCACLDKQRPCQSYN GKQVYKETCISFIKD FLAKNLIVIGIAFGLAVIEILGLVFSMVLYCQIGNK |
| Amino acids | 1 to 237 |
| Tags | GST tag N-Terminus |

Specifications

Our [Abpromise guarantee](#) covers the use of **ab159720** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|-------------------------|---|
| Applications | Western blot ELISA |
| Form | Liquid |
| Additional notes | Protein concentration is above or equal to 0.05 mg/ml. Protein previously labeled as TM4SF3. |

Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.


pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

| | |
|------------------------------|--|
| Tissue specificity | Gastric, colon, rectal, and pancreatic carcinomas. |
| Sequence similarities | Belongs to the tetraspanin (TM4SF) family. |
| Cellular localization | Membrane. |

Images



ab159720 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human Tspan-8 protein
(ab159720)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors