

Product datasheet

Recombinant Human TXNDC15 protein ab194018

Description

Product name	Recombinant Human TXNDC15 protein	
Purity	> 95 % SDS-PAGE. Purity is greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.	
Endotoxin level	< 1.000 Eu/μg	
Expression system	HEK 293 cells	
Accession	Q96J42	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	VEVAEESGRLWSEEQPAHPLQVGAVYLGEEELLHDP MGQDRAAEEANAVL GLDTQQGDHVMVLSVIPGEAEDKVSSEPSGVTCGAGG AEDSRCNVRESLFS LDGAGAHFPDREEEYYTEPEVAESDAAPTEDSNNTES LKSPKVNCEERNI TGLENFTLKILNMSQDLMDFLNPNNGSDCTLVLFYTPWC RFSASLAPHFNS LPRAFPALHFLALDASQHSSLSTRFGTVAVPNILLFQG AKPMARFNHTDR TLETLKIFIFNQTGIEAKKNVVVTQADQIGPLPSTLIKSVD HHHHHH	
Predicted molecular weight	33 kDa including tags	
Amino acids	33 to 321	
Tags	His tag C-Terminus	
Additional sequence information	Extracellular domain.	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab194018** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

	HPLC
Form	Lyophilised

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. pH: 7.40 Constituents: 0.87% Sodium chloride, 99% Phosphate Buffer
Reconstitution	Reconstituted protein solution can be stored at 4-7 deg.C for 2-7 days. Aliquots of reconstituted samples are stable at < -20 deg.C for 3 months.

General Info

Sequence similarities	Contains 1 thioredoxin domain.
Cellular localization	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors