

Product datasheet

Recombinant Human WDSUB1 protein ab165627

1 Image

Description

<b>Product name</b>	Recombinant Human WDSUB1 protein
<b>Expression system</b>	Wheat germ
<b>Protein length</b>	Full length protein
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Sequence</b>	<p>MVKLIHTLADHGDDVNCCAFSFSLLATCSLTKTIRLYSL  RDFTELPHSPL  KFHTYAVHCCCFSPSGHILASCSTDGTTVLWNTENGQ  MLAVMEQPSGSPV  RVCQFSPDSTCLASGAADGTVVLWNAQSYKLYRCGS  VKDGSLAACAFSPN  GSFFVTGSSCGDLTVWDDKMRCLHSEKAHDLGITCC  DFSSQPVSDGEQGL  QFFRLASCGQDCQVKIWSFTDILARRTEHQLKQFTE  DWSEEDVSTWLC  AQDLKDLVGIFKMNNIDGKELLNLTKESLADDLKIESLG  LRSKVLRKIEE  LRTKVKSLSSGIPDEFICPITRELMKDPVIASDGYSYEK  EAMENWISKKK  RTSPMTNLVLPsAVLTPNRTLKMAINRWLETHQK</p>
<b>Amino acids</b>	1 to 384
<b>Tags</b>	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab165627** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	Western blot ELISA
<b>Form</b>	Liquid
<b>Additional notes</b>	Protein concentration is above or equal to 0.05 mg/ml.

## Preparation and Storage

---

### Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

## General Info


---

### Relevance

WD repeat-containing proteins are those that contain 4 or more copies of the WD-repeat (tryptophan-aspartate repeat), a sequence motif approximately 31 amino acids long, that encodes a structural repeat. There are two isoforms of WDSUB1.

## Images

---



ab165627 on a 12.5% SDS-PAGE stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors